

Analysis of broad occupancy patterns

Epigenomics Data Analysis Workshop

Stockholm, 24 November 2020

Agata Smialowska



NBIS, SciLifeLab, Stockholm University




Point-source vs. broad peak detection

Sequence-specific binding (TFs)

Distributed binding (histones, RNAPol2)



How to define sufficient read depth?



detection sensitivity - recovery of true enrichment region




percent increase in enriched regions recaptured when an additional 1 million reads are sequenced



'sufficient depth' - the sequencing depth at which the percent gain per 1 million additional sequence reads falls below 1%

Effect of sequencing depth on site recovery




Sequencing depth

Chromatin
Remodellers
Histone marks



mixed signal

Chromatin
Remodellers
Histone marks
RNA polymerase II



broad signal

Human:
3.2Gb

H3K36me3: 25 M

H3K27me3: 35 M

H3K9me3: >55 M

Drosophila:
300Mb

H3K36me3: 11 M

H3K27me3: 20M

H3K9me3: 20 M

No clear guidelines for mixed and broad type of peaks

Source: The ENCODE consortium; Jung et al, NAR 2014

Peak calling

Chromatin
Remodellers
Histone marks



mixed signal

Chromatin
Remodellers
Histone marks
RNA polymerase II




broad signal

- **MACS2** in broad mode: window-based method using local statistics
- **csaw**: Scoring in moving windows
- **Epic2** (SICER): tendency of histone modifications to cluster to form the domains. This method identifies islands as *clusters* of enriched windows. Islands, rather than individual windows of fixed length, are the fundamental units of interest


Effect of sequencing depth on regions detected by various algorithms

b



Percent of recaptured enriched regions
All enriched regions



Top 20% enriched regions




Comparison of enriched regions detected by various algorithms



55M
human


Comparison of enriched regions detected by various algorithms



Exercise

- MACS2 in `broad peak` mode;
- csaw: Detection of differentially bound regions in ChIP-seq data with sliding windows, with methods for normalization and proper FDR control;
- Low sequencing depth data; H3K79me2 (transcribed regions of active genes)

MACS2 vs. csaw



Resources for broad region analysis

- <https://omictools.com/peak-calling-category>
- <https://www.encodeproject.org/chip-seq/histone/>