

Epigenomics methods at NCI



Orlando Contreras-López, Ph.D.

Leader Genomic Applications Lab Team
orlando.contreras@scilifelab.se

Overview

The **National Genomics Infrastructure (NGI)**

DNA methylation

Omni-C and Micro-C

Omni-ATAC

Nanopore Sequencing

NGI Sweden Nodes

 NGI stockholm



The **SNP&SEQ**
Technology Platform




Uppsala Genome
Center - **UGC**



<http://snpseq.medsci.uu.se/>

<http://ugc.igp.uu.se/>

NGI Stockholm




<https://ngisweden.scilifelab.se/resources/ngi-stockholm-status/>


NGI Workflow




Meetings



QC & Library
Preparation




Sequencing



Data analysis
& delivery

We provide **guidelines and support** for sample collection, study design, protocol selection and bioinformatics analysis

NGI Workflow



NGI Workflow



Meetings

- 1) Here are a few challenges that users often face:
 - Not knowing enough about their own sample
 - Lacking the necessary tools for quality control
 - Needing more bioinformatics capabilities
 - Feeling lost in communication/literature

- 2) NGI's tools
 - Knowledge gathered from diverse sample types
 - Providing guidance and sample QC service, clear thresholds and requirements.
 - Offering basic data QC and connecting to NBIS
 - Accessing the website and scheduling one-to-one meetings



De novo sequencing



DNA extractions

DNA resequencing



Environmental sequencing

Epigenetics



Genotyping

Single-cell



Spatial

Transcriptomics



User-prepared libraries

Proteomics with NGS readout



Overview

SciLifeLab

The National Genomics Infrastructure

DNA methylation

Omni-C and Micro-C

Omni-ATAC

Nanopore Sequencing

DNA methylation analysis

Library preparation for whole genome methylation sequencing

1. NEBNext® Enzymatic Methyl-seq (EM-Seq)
2. SPLinted Ligation Adapter Tagging (SPLAT) for Whole Genome Bisulphite sequencing

*Developed at NGI!**

1. Illumina EPIC Arrays.

*Amanda Raine, Erika Manlig, Per Wahlberg, Ann-Christine Syvänen, Jessica Nordlund, SPLinted Ligation Adapter Tagging (SPLAT), a novel library preparation method for whole genome bisulphite sequencing, *Nucleic Acids Research*, Volume 45, Issue 6, 7 April 2017, Page e36



Overview

SciLifeLab

The National Genomics Infrastructure

DNA methylation

Omni-C and Micro-C

Omni-ATAC


Nanopore Sequencing

Epigenetics - HiC


Chromosome territories




A/B compartments



TADs



Loops




Resolution
(Seq Depth)




doi.org/10.1016/j.tibs.2019.03.001

HiC: Restriction Enzyme




Dovetail Genomics

Epigenetics - HiC




doi.org/10.1016/j.tibs.2019.03.001

Epigenetics - HiC



Resolution
(Seq Depth)


HiC: Restriction Enzyme



Lysis



No Fragmentation

DNase/MNase digestion



Dovetail Genomics

Omni-C vs Micro-C



- **Enzymatic titration**
- **Sample quality**
- **Frozen cells**
- **Frozen nuclei**
- **Genome Assembly**
- **Frozen Animal tissues**
- **Frozen Plant tissues/nuclei**

10.1016/j.cell.2013.09.011

-C steps at NGI

→ Contact and Planning

- ◆ Study design
- ◆ Aims and Biological Question.
- ◆ Expectations and Limitations
- ◆ Turnaround Times
- ◆ Paperwork and Sample delivery



→ 200M reads / 1Gb genome to start with

→ Around 6500 kr / Sample

Overview

SciLifeLab

The National Genomics Infrastructure


DNA methylation

Omni-C and Micro-C

Omni-ATAC

Nanopore Sequencing


Epigenetics - ATAC-seq



- Nucleosome mapping
- Transcription factor binding analysis
- Novel enhancer identification
- Evolutionary studies
- Comparative epigenomics

Illumina


Epigenetics - Omni ATAC




[corceslab.com/pages/Resources/
10.21769/BioProtoc.4294s/](https://corceslab.com/pages/Resources/10.21769/BioProtoc.4294s/)

Epigenetics - Omni ATAC

FastQC: Sequence Counts




Preseq: Complexity curve




Omni ATAC has a higher number of unique reads, compared to regular ATACseq

Epigenetics - Omni ATAC



Omni ATAC generates fewer reads that map to mitochondrial DNA, compared to regular ATACseq



TSS enrichment between comparable for all samples

All annotated genes have been normalized to the same size.

- Green: -3.0Kb upstream of gene to TSS
- Yellow: TSS to TES
- Pink: TES to 3.0Kb downstream of gene

Omni ATAC at NGI

→ Contact and Planning

- ◆ Study design
- ◆ Aims and Biological Question.
- ◆ Expectations and Limitations
- ◆ Turnaround Times
- ◆ Paperwork and Sample delivery



→ Sequence > 50 M PE reads

→ Around 2000 kr / Sample

NovaSeq 6000

Overview

SciLifeLab

The National Genomics Infrastructure


DNA methylation

Omni-C and Micro-C



Omni-ATAC

Nanopore Sequencing

Oxford Nanopore Technologies (ONT)




- A protein channel (nanopore) anchored to a non-conductive membrane
- A voltage can be applied across the membrane
- This causes ions to travel through the protein channel, creating an electrical current



ONT new possibilities

- **Sequencing native DNA or RNA** avoids bias introduced by amplification steps
- **RNA/DNA Methylation calling**
- Long reads advantages:
 - **haplotyping**
 - analysis of repetitive regions
 - phasing
 - **structural variants**
- Very simple library preps and easy instrument operation

ONT Adaptive Sampling



- ~150X coverage of the ROI (2Mb)
- 20-40X coverage for the rest of the genome
- No special Lib Prep
- Dependent on .bed file design

<https://ngisweden.scilifelab.se/2024/02/tech-note-adaptive-sampling-targeted-oxford-nanopore-long-read-sequencing/>

ONT coming soon

Future methods:

- **Direct RNA-seq**
- **Single Cell 10X Genomics**
 - 3' Gene Expression
- **Spatial Transcriptomics 10X Genomics**
 - Visium Spatial
- **Pore-C: HiC + ONT**



→ **Around 2300 kr / Sample**



N50 (Kb)	0.41	32.4
Output (Gb)	44.8	34.8
Reads (M)	111.2	3.6

Takeaway

- Make use of the SciLifeLab platforms!
 - National Genomics Infrastructure
 - <https://ngisweden.scilifelab.se/>
 - National Bioinformatics Infrastructure
 - <http://nbis.se/>
- Pick a protocol that matches your biological question
- Think carefully about your project plan
- Ensure that you understand your data!
- **Ask us!** -even if we don't offer the service, we might...

